

HFP-SS90/160

Mechanical double arm rotating pendant anesthesia pendant

Describe

Applicable areas: operating room, emergency room, intensive care unit, recovery room. The installation adopts the ceiling hanging mode, which can rotate within 340 ° range, and accurately adjust the position of the terminal box. The height of the equipment enables the medical staff to raise their hands and be handy. The suspension tower cantilever is made of imported high-strength aluminum alloy by one-step extrusion, and the surface is treated with primary oxidation and rust prevention. Brake: optional friction brake or pneumatic brake.

Standard Configuration

- * 1 oxygen (O2)
- * 1 air
- * 1 nitrous oxide (N2O)
- * 1 vacuum suction (VAC)
- * 1 jet exhaust emission (AGSS)
- * 6 x 10A power outlets
- * 1 network socket
- * 1 grounding terminal

Annex

- * 2 pallets
- * 1 drawer
- * 1 infusion stand



SHANGHAI HUIFENG MEDICAL INSTRUMENT CO.,LTD.

Company Address: No. 885, Qixing Road, Nичeng Town, Pudong

New Area, Shanghai

Zip code: 201306

Sales call: 021-58082099, 58085105, 58084992, 58086595

Sales Fax: 021-58085105

After-sales service telephone: 021-58084992

After-sales service fax: 021-58085105

Email address: lisa@huifeng-co.com

Website: www.sh-huifeng.com

The company has passed ISO 9001, ISO 13485 quality management system certification